


# Inserting the Triple Flange or Quad Flange Earplugs

Losing your hearing is not “just part of the job” wearing your hearing protection is!



Place the stem of the TF earplug into the seating device or firmly grasp the stem of the QF earplug.



Gently push or twist earplug into ear canal.



The last flange of the TF earplug should be flush with the ear canal opening. The last flange of the QF earplug may or may not be flush with the ear canal opening. Tension should be felt when lightly pulling on the earplug stem.

**Conserve fighting strength**  
**Conserve your hearing**



## WEAR

- Own voice will sound low toned. Sounds will be muffled
- Twist rather than pull earplugs out of ear
- Re-insert periodically if worn for long periods

## CARE

- Wash earplugs with warm soapy water
- Allow earplugs to completely dry before placing into the ear or earplug case
- Store earplugs in an earplug case to minimize loss

## ORDERING

### Triple Flange

- NSN Small 6515-00-442-4821
- NSN Large 6515-00-467-0092
- NSN Medium 6515-00-442-4818

### Quad Flange

- NSN 6515-01-492-0443

### Earplug Case

- NSN earplug case 6515-01-100-1674



For more information on noise and hearing loss visit  
<http://chppm-www.apgea.army.mil/hcp/>